

Fango Stirrers RW / NFR

Medical Device Class I



Fango Stirrer RW 44, mounted on Heating Cabinet WS 14-6044 (Mounting example)

Intended Use

Type RW: Heating and tempering of fango paraffin in physiotherapy

Type NFR: Heating and tempering of natural fango in physiotherapy



Equipment

- Self-supporting housing made of stainless steel with structured surface
- Aluminium or stainless-steel vessel (Type NFR: stainless steel vessel)
- Top lid with red, green or black handle bar, lid fastened with two hinges
- Energy saving insulation
- Electronic temperature control with digital temperature display
- Safety thermometer, displays the temperature even when the device is switched off
- Indicator lamps for heating, stirring and high temperature treatment
- Continuous stirring, interval stirring
- Heated ball valve 1 ¼" for fango paraffin (Type RW)
- Ball valve 2" for natural fango (Type NFR)
- Motor protection relay

Technical Data

Device Article Number	Filling Capacity (kg)	Heating Capacity (W)	Nominal Voltage (V)	Outer Dimensions (mm) L x T x H	Temperature Range
RW 44 E 0020.100 AI	40	1000	230	540 x 540 x 590	30 °C - 110 °C 135 °C
RW 64 E 0022.100 AI 0022.200 VA	60	1000 (AI) 2000 (VA)	400	620 x 620 x 655	30 °C - 110 °C 135 °C
RW 84 E 0024.100 AI 0024.200 VA	80	2000	400	680 x 680 x 682	30 °C - 110 °C 135 °C
RW 104 E 0025.100 AI 0025.200 VA	100	3000	400	720 x 720 x 660	30 °C - 110 °C 135 °C
NFR 64 E 0022.300	60	1000	400	620 x 620 x 655	Room temperature to 65 °C
NFR 84 E 0024.300	80	2000	400	680 x 680 x 682	Room temperature to 65 °C
NFR 104 0025.300	100	3000	400	720 x 720 x 660	Room temperature to 65 °C

All dimensions are given in SI units